

**BRIKEN Commissioning November 2016**  
**Running Gasific-7.0 acquisition program**  
(07/11/2016)

**Computer:** IFICdaq.riken.jp (IP: 10.32.6.166)

**User:** dacq

**Password:** BRIKEN2016

**Working directory:** /data/1611Commissioning/161107Parasitic/

**Configuration file:** 161107Conf\_BrikenFull.xlsx

### 1) Start Gasific-7.0

- Make sure that you are on **Desktop 1**
- You should see four terminals labeled *AcqRead*, *GUI*, *RootDisplay*, and *ConfClock*. **NEVER KILL THESE TERMINALS.**
- Check that you are in the right working directory on the 4 terminals
- In the *ConfClock* window type **ConfClock\_external** and press enter to start the common clock for synchronization with AIDA and BigRIPS. The window becomes **GREEN**. Unless the VME crates are reset there is no need to repeat the command even if you stop Gasific-7.0. If the *ConfClock* terminal is yellow you will not be able to synchronize.
- Type in the corresponding terminal the commands **AcqRead**, **GUI\_DAQ**, and **RootDisplay**, and press enter in order to launch the Gasific-7.0 processes. The **MainWindow** panel and the **ROOT canvas** will open.
- To stop Gasific-7.0 processes: ctrl+C in all of the three terminals *AcqRead*, *GUI*, and *RootDisplay*. If the data acquisition stacks or you don't have control on the *GUI* or *ROOT display* you will need to do that.

### 2) Load the proper configuration file

In the **MainWindow** panel:

- Click on **Config Module** tab
- Click on **Load Config** button
- A menu with possible configuration files will appear. Select the right configuration file (\*.xlsx), and click the **Open** button
- Wait until the file is read (monitor the progress in the *GUI* terminal) and click the **Send to Hard** button. **Be patient** it takes about 1 minute to load all parameters in the digitizers (monitor the progress in the *AcqRead* terminal, wait until the line "List of controllers defined" followed by 4 additional lines appear)

### 3) Start the synchronization with AIDA/BigRIPS

Needed in principle only once. Needed if any of the three DACQs crashes (start at point 1). Not needed for running standalone.

In the **DAQ Control** tab:

- Click **Start DAQ**
- Click the **Start BRIKEN Sync** button. Enter in **Server host** the address **10.32.0.12** and click **OK**

- Click **Send Sync pulse** to reset all timestamps to 0.
- Click **Stop DAQ**

#### 4) Start/stop acquisition

In the **DAQ Control** tab

- You may want to click the **Reset** button to clear the *DAQ Statistics* counters

**Without list-mode data:**

**Start:**

- Click on **Start Online** button
- Click on **Start DAQ** button

At this point you should see that the *DAQ Statistics* counters are being updated.

**Stop:**

- Click on **Stop DAQ** button
- Click on **Stop Online** button

**With list-mode data:**

**Start:**

- Click on the **Start File** button. You are prompted for a **File name**. Use the format **YYMMDD\_HHMM\_\*** (the extension **.dlt** will be added by the system). **DO NOT ENTER A RUN NUMBER**. The *Run number* is automatically incremented by the system. Instead check that the *Run number* has been effectively incremented respect to the last written file (if this is not the case then provide a run number higher than the last one). Click the **OK** button.
- Click on the **Start Online** button
- Click on the **Start DAQ** button

**Stop:**

- Click on the **Stop File** button in order to close the current **dlt** file
- Click on **Stop DAQ** button
- Click on **Stop Online** button

#### 4) Visualize online spectra

In the **Online** tab:

- Click on **Raw, Calibrated, Correlated, or Groups** sub-tabs to see the list of corresponding types of spectra to display
- **Click on the name** of the spectra that you want to visualize. It will be displayed in the *RootDisplay* window (*Canvas*)
- On the *Canvas* you can select **Auto Update, Log Scale or Overlay**. You can also use the usual interactive ROOT commands (zoom, change line color, fit, ...).

#### 5) Save and clear online spectra

In the **Online** tab:

- Click on the **General** sub-tab
- Click on **Save online histograms** to save all online spectra into a root file. You are prompted for a file name. Enter it using the format **YYMMDD\_HHMM\_\*.root** and click **Save**. **DO**

**NOT USE ENTER.** Please **verify** that the file has been effectively written on disk.

- Click on **Clear online histograms** to delete all online spectra. You are asked to confirm the operation. **THINK BEFORE YOU ERASE THE HISTOGRAMS.**

**Troubleshooting: Call an expert!** (Jorge, Alvaro, Jose)

### LIST OF HISTOGRAMS:

**Raw:** V1A1C1, ..., V1A7C16, V2A1C1, ..., V2A7C8

**Groups:** GV1A1, ..., GV1A7, GV2A1, ..., GV2A7

**Calibrated:**

He001, ..., He140, D4Black, D4Red, D4Green, D4Blue, G7Black, G7Red, G7Green, G7Blue, F11\_PL\_R, F11\_PL\_L, Si\_T, Si\_B, V\_PL\_T, V\_PL\_B, AIDA\_PL, V\_NaI\_T, V\_NaI\_B, V\_NaI\_R, V\_NaI\_L, Pulser

**Correlated:**

BRIKEN, Ring1, ..., Ring7, D4, G7, Clovers  
RateBRIKEN, RateR1, ..., RateR7, RateD4, RateG7, RateG7Black, RateG7Green, RateClovers  
RateF11R, RateVPLT, RateVPLB, RateAIDAPL, RateSiT, RateSiB  
Corr\_Br\_BR, Corr\_R1\_R2, Corr\_BR\_D4, Corr\_BR\_G7, Corr\_BR\_Clovers, Corr\_D4\_G7,  
CorrG7BlackGreen, Corr\_BR\_F11R, Corr\_BR\_VPLT, Corr\_BR\_VPLB, Corr\_BR\_AIDAPL,  
Corr\_BR\_SiT, Corr\_BR\_SiB, Corr\_F11R\_Clovers, Corr\_F11R\_VPLT, Corr\_F11\_VPLB

## RUN PROTOCOL:

- 1) Start a new measurement:
  - a. **Start File** (enter dlt filename: **YYMMDD\_HHMM\_\***)
  - b. **Start Online**
  - c. **Start DAQ**
  - d. Reset DAQ statics counters
- 2) Check that the DAQ statistics counters are counting. Check that the BRIKEN spectrum is incrementing. Check that the new \*.dlt file has been opened and is incrementing.
- 3) Make the corresponding entry in the ELOG: Start time, BRIKEN DLT filename, BRIKEN run number, AIDA filename, AIDA run number, Run info (Beam, BigRIPS setting, ...)
- 4) Fill the run-sheet
- 5) Check the important histograms: BRIKEN, Clovers, F11\_PL\_R, Si\_T, V\_PL\_T, AIDA\_PL, .... Check the Rate\* histograms. Check the Corr\* histograms.
- 6) Stop the measurement:
  - a. **Stop File**
  - b. **Stop Online**
  - c. **Stop DAQ**
  - d. **Save online histograms** (enter ROOT filename: **YYMMDD\_HHMM\_\*.root**)
- 7) Check that the \*.root file has been written.
- 8) Fill the information in the ELOG (Stop time, BRIKEN run number, AIDA run number, Rates, additional info)
- 9) Complete the run-sheet info
- 10) **Clear online histograms**